



REHASENSE

LIGHTWEIGHT  
FUNCTIONALITY



# SERVER

LIGHTWEIGHT ROLLATOR

The stylish Server rollator has a unique triangle profile ensuring a steady frame design. It is among the lightest and easiest to turn rollators on the market in the 150 kg user category.

# SERVER

LIGHTWEIGHT ROLLATOR



REHASENSE



## FOLDABLE

The Server folds easily and stays that way with a simple locking system that doubles up as an ergonomically shaped carrying handle.

## STEADY AND DURABLE


The frame and seat are tested and approved for maximum user weight of 150 kg. It has light and reliable brakes with TPE covered unbreakable wheels. It has a practical removable textile shopping bag.

## HEIGHT ADJUSTABLE

The Server is available in 3 sizes – small, medium and large – with a range of settings within these sizes. The push handles are easy height adjustable in a range of 60 - 102 cm.



Smooth adjustable handle height and lockable handbrakes.



Seat heights 50, 55 and 62 cm.



Maintenance free brakes and TPE-covered unbreakable wheels.



Unique triangle aluminium profile ensures the durability up to 150 kg.



Stays firmly closed with the lock and is easy to open with the push button.



Practical easy-to-remove shopping bag with special designed easy-to-go clips.

## TECHNICAL INFORMATION

Article number (example Wine Red)  
SRWRL600 (Wine Red Large)  
SRWRM550 (Wine Red Medium)  
SRWRS500 (Wine Red Small)

### Colours Small / Medium / Large:

WR Wine Red  
BO Mocca Brown  
CH Champagne

### Colours only Large:

TU Turquoise  
LIL Lila  
SG Sober Grey  
BK Phantom Black  
WH White

### Seat Height:

Small 50 cm  
Medium 55 cm  
Large 62 cm

### Seat width: 46 cm

### Measurements: (width x depth x height)

Small 60 x 66 x 60–72 cm  
(users height 125–160 cm)  
Medium 60 x 66 x 66–86 cm  
(users height 135–170 cm)  
Large 60 x 66 x 74–102 cm  
(users height 150–200 cm)

### Transportation measurements folded:

(depth x height x width)  
Small: 66 x 72 x 22,5 cm  
Medium: 66 x 74 x 22,5 cm  
Large: 66 x 79 x 22,5 cm

### Shopping bag: (width x depth x height)

32 x 15 x 28 cm

### Weight:

Small 6,7 kg (+ basket 370 g)  
Medium 6,8 kg (+ basket 370 g)  
Large 7,0 kg (+ basket 340 g)

Max. user weight: 150 kg

Frame material: aluminium

## COLOURS

SRBOL600 Mocca Brown

SRWHL600 White

SRTUL600 Turquoise

SRCHL600 Champagne

SRBKL600 Phantom Black

SRLIL600 Lila

SRWRL600 Wine Red

SRSGL600 Sober Grey

Server is CE-marked and fulfills the requirements of EU directive 93/42 EEC for medical devices.

Rehasense sp. z o.o.  
Sulejowska 45  
97-300 Piotrków Trybunalski  
Poland

Rehasense Benelux  
Markeweg 81  
8398 GM Blesdijke  
Netherlands

Rehasense Deutschland  
Briener Strasse 25  
47533 Kleve  
Deutschland

Copyright © Rehasense® 2017