



REHASENSE

CARBON FIBRE
ELEGANCE AND
PERFORMANCE



ATHLON SL

CARBON FIBRE TECHNOLOGY

Athlon SL – carbon fibre composite rollators highlight the elegance and performance of the Rehasense range of rollators. Lightweight from 5 kg, **Athlon SL** can be used for users up to 150 kg.

ATHLON SL

CARBON FIBRE TECHNOLOGY



REHASENSE



FOLDABLE AND EASY TO CARRY

Athlon SL is easy to fold and stays folded with an easy-to-use locking system. Ergonomically designed handle grips make Athlon easy to transport, if needed.

STEADY AND DURABLE


The frame and seat are tested and approved for the maximum user weight of 150 kg. Light and reliable brakes with TPE-covered unbreakable wheels, create a steady and safe driving feeling.

HEIGHT ADJUSTABLE

Athlon SL push handles height adjustable to 10 positions, from 63 cm to 102 cm. User comfort is ensured with three different seat heights - 62, 55 and 50 cm - as well as optimal frame design and ergonomic push handles.



Smooth adjustable handle height and lockable handbrakes.



Seat heights 62, 55 and 50 cm.



Maintenance free brakes and TPE-covered unbreakable wheels.



Unique triangle carbon fibre profile ensures the durability up to 150 kg.



Stays firmly closed with the lock and is easy to open with the push button.



Practical easy-to-remove shopping bag with special designed easy-to-go clips.



Ergonomically designed handle grips make **Athlon SL** easy to transport, if needed.



Athlon SL is easy to fold and stays folded with an easy-to-use locking system.



All sizes can be fitted with soft comfort wheels.

ACCESSORIES

Athlon SL can be fitted with a broad range of accessories from Rehasense.

TECHNICAL INFORMATION

Colors:



Max. user weight: 150 kg

Item no:

Large
CRBKL600 (black)
CRWRL600 (red)
CRSGL600 (anthracite)
CRWHL600 (white)

Medium
CRBKM550 (black)
CRWRM550 (red)

Small
CRBKS500 (black)

Measurements:

(width x depth x height)

Large: 61 x 68 x 74-102 cm

Medium: 61 x 68 x 66-86 cm

Small: 61 x 68 x 63-79 cm

Seat width: 46 cm

Transportation
measurements folded:

(depth x height x width)

Large: 68 x 79 x 22,5 cm

Medium: 68 x 74 x 22,5 cm

Small: 68 x 72 x 22,5 cm

Weight:

Large: 5,2 kg (+ bag 370 g)

Medium: 5,1 kg (+ bag 370 g)

Small: 5,0 kg (+ bag 340 g)

Shopping bag:

(width x depth x height)

32 x 15 x 28 cm

Max. load 5 kg.

Athlon SL is CE-marked and fulfils the requirements of EU directive 93/42 EEC for medical devices.

Copyright © Rehasense® 2018

Rehasense sp. z o.o.
Sulejowska 45g
97-300 Piotrków Trybunalski
Poland

Rehasense Benelux
Markeweg 81
8398 GM Blesdijke
Netherlands

Rehasense Deutschland
Briener Strasse 25
47533 Kleve
Deutschland